





# Compex®

www.compex.info

# Compex®

- Position confortable
- Bequeme Haltung
- Posizione confortevole
- Comfortable position
- Posición confortable
- Comfortabele positie

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## INTRODUCTION 121

## WARM-UP CATEGORY 122

Preparation for seasonal endurance activity (eg: hiking, bicycle touring)	123
Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, mountain biking, mountain stage in cycling, cycling pursuit on track, 200 metres swimming, slalom)	123
Prevention of contractures in the muscles of the back of the neck for cyclists	123
Muscular and circulation conditioning before physical activity	123
Use of Potentiation program to optimise the effects of explosive strength (sprints, jumps, squash, football, basketball, etc.) immediately before specific training or competition	124

## SPORT CATEGORY 125

Preparation for a cyclist training three times a week who wishes to improve endurance performance	126
Preparation for a runner training three times a week who wishes to progress in endurance (half-marathon, marathon)	126
Pre-season preparation of lactic capacity for a resistance sport with three active training sessions per week (800 metres, track cycling, etc.)	127
Preparation for a cyclist training three times a week who wishes to improve his power	127
Preparation for a swimmer training three times a week who wishes to improve his swimming power	128
Pre-season preparation for a team sport (football, rugby, handball, volleyball, etc.)	128
Maintaining results achieved during preparation for team sports during the competitive period (football, rugby, handball, volleyball, etc.)	129
Pre-season preparation of the explosive strength of the quadriceps for an athlete training three times a week (long jump or high jump, sprinting, etc.)	129
Preparation for a thrower (javelin, shot putt) or tennis player (service, smash) who wishes to increase explosive strength through plyometric training	130
Muscle volume for a body-builder	130

Muscle volume for a body-builder training three times a week	131
Muscle volume for a body-builder training five times a week	131

## CROSS-TRAINING CATEGORY 132

Resuming activity after a break in training: variable muscle preparation	133
Preparation for a long-distance runner who also wants to improve his speed at the end of the race	133
Improving the user's capacity for intense efforts without increasing the risks of injury	134
Competition preparation for an amateur practising a team sport	134

## CONDITIONING CATEGORY 136

Abdominal muscle workout for golfers	137
Strengthening the lower back muscles in rowers	137
Improving stabilisation of the upper body muscles in triathletes	137
Preventive stimulation of the lateral peroneal muscles to avoid or reduce the risk of recurrence of ankle sprain in footballers	138

## RECOVERY CATEGORY 139

Practical rules of use for the Undulated massage 1 and Undulated massage 2 programs	141
Increasing speed of recovery after intense competitive events, such as marathons	142
Treatment of uncomfortable muscular tensions in the back of the neck	142
Treatment of a localised sensation of heaviness or an occasional state of fatigue	142
Returning to a state of well-being after stressful situations	143
To eliminate unpleasant body sensations and recover an optimal state of relaxation	143
To restore the feeling of muscle equilibrium and soothe undesirable manifestations	143
Use of the Regeneration program to eliminate muscle fatigue more quickly (cross-country running, football, basketball, tennis, etc.) and restore good muscular sensations more quickly	143
Prevention of the feeling of heaviness in the legs	144

Preparation for a footballer wishing to optimise the effects of active stretching of the hamstrings muscles	144
---	-----

**PAIN RELIEF CATEGORY 146**

Neuralgia of the upper limb (brachial neuralgia)	148
Chronic muscular pain (polymyalgia)	148
Contracture (eg: localized contracture in external side of the calf)	148
Chronic muscular pain in the back of the neck (cervical pain)	148
Muscular pain in the thoracic region (thoracic back pain)	149
Muscular pain in the low back region (low back pain)	149
Sharp and recent muscular pains affecting a muscle in the low back (lumbago)	149
Chronic elbow pain (epicondylitis = tennis elbow)	149
Torticollis (recent and acute back of the neck muscle pains)	150

**REHABILITATION CATEGORY 151**

Atrophy of the quadriceps as a result of trauma	152
Development of the latissimus dorsi to treat and prevent tendinous shoulder pain (rotator cuff syndromes)	152
Development of the abdominal belt to prevent pain in the lumbar region (low back pain)	152
Development of the low back muscles to prevent pain in the lumbar region (low back pain)	152
Strengthening of the lateral peroneus muscles after ankle sprain	153

**TEST CATEGORY 154**

Locating a motor point (eg: locating the motor points of the vastus internus and the vastus externus of the quadriceps)	154
---	-----

## INTRODUCTION

The applications that follow are given by way of example. They provide a better understanding of how electrostimulation sessions can be combined with voluntary activity. These protocols will help you to determine the best procedure to follow depending on your needs (choice of the program, muscular group, duration, placement of the electrodes, body position).


You can of course choose a muscular group different from the one indicated, depending on the activities you practice or the body area you wish to stimulate.

The specific applications provide you information on the placement of electrodes and the stimulation position to adopt. This information is represented by numbers (placement of electrodes) and letters (stimulation position). They refer to the pictures for the placements of electrodes and the pictograms for the body positions appearing on the flap of the manual's cover.

When the Compex muscle stimulator **is not equipped** with the **Mi-SENSOR**, you must select in the device the muscular group indicated on the figurine. When no figurine appears (–), it is not necessary to choose the muscular group: the selection is automatic.

When the Compex muscle stimulator **is equipped** with the **Mi-SENSOR**, selection of the muscle group is automatic. Consequently, it is unnecessary to refer to the small pictogram.

Example :

For the four basic training programs, **Endurance, Resistance, Strength and Explosive strength**, you are strongly advised to consult the Training planner available in the DVD. An interactive questions and answers system gives you access to a personalized training plan.

# Warm-up Category

Programs	Effects	Uses	Stimulation energies	mi functions
<b>Capillarization</b>	Very marked increase in blood flow Development of capillaries	In the period before competition, for those participating in endurance or resistance sports  To improve endurance in those whose fitness level is not high	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Toning massage</b>	Activation of blood circulation Recovery of muscular contractile properties Invigorating effect	To prepare the muscles in an ideal manner before an unusual/one-time physical activity	Progressively increase the stimulation energies until pronounced muscle twitches are obtained  Make sure that the stimulation energy is sufficient in order to impose significant muscle contractions	mi-SCAN mi-RANGE
<b>Potentiation</b>	Increased speed of contraction and power gain Less nervous effort to attain maximum strength	For optimum muscle preparation immediately before competition	Maximum tolerable energy (0-999)  The energy is adjusted during the phases of the program when the muscle relaxes and reacts with simple muscle twitches	mi-SCAN

## Warm-up

### Specific applications

#### Preparation for seasonal endurance activity (eg: hiking, bicycle touring)

Sporting physical activities that last a long time (walking, cycling, cross-country skiing, etc.) require strong muscles with a good capillary circulation so that the muscle fibres can be properly oxygenated. When a long-duration physical activity is not practised, or not practised regularly enough, the muscles lose their capacity to effectively consume oxygen and the capillary system becomes rarefied. This lack of muscular quality makes exercising uncomfortable and limits the speed of recuperation, and is responsible for many unpleasant side effects such as drowsiness, contractures, cramp and swelling.

To restore to muscles their endurance and ensure full activation of the capillary system, the Compex offers a very efficient mode of stimulation.

Cycle duration: 6-8 weeks, before starting the hike, 3 x/week, alternating muscular groups

Program: **Capillarization 8**  and **Capillarization 25** 

#### Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, mountain biking, mountain stage in cycling, cycling pursuit on track, 200 metres swimming, slalom)

Sports that require maximum effort lasting between 30 seconds and 5 minutes are the so-called resistance sports. Rapid fibres must work at a capacity close to their maximum and be capable of maintaining this intensive work throughout the event without weakening, in other words the rapid fibres must be resistant.

The **Capillarization** program, which produces a very significant increase in the blood flow in the muscles, brings about development of the intramuscular capillary system (capillarization). This growth of the capillaries works best if it occurs around the rapid fibres. In this way, the latter's exchange surface with the blood increases, allowing an improved supply of glucose, better diffusion of oxygen and quicker evacuation of lactic acid. Capillarization therefore allows the rapid fibres to be more resistant and to maintain their optimum capacity over a longer period.

However, prolonged or too frequent use of this program may induce a modification of rapid fibres into slow fibres, thereby risking a decline in performance for strength and speed sports. It is therefore important to follow closely the recommendations below concerning the specific application in order to benefit from the positive effects of this treatment.

Cycle duration: 1 week before competition, 2 x/day, with a 10 minute break between the 2 sessions

Program: **Capillarization 8** 

#### Prevention of contractures in the muscles of the back of the neck for cyclists

During voluntary training or competitions, some athletes experience the problem of contractures in support muscles (for example, the back of the neck for cyclists). Use of the **Capillarization** program reduces, and may resolve, this problem. The marked increase in blood circulation and development of the capillary system improve the oxygenation of fibres and their exchanges with the blood.

The muscle therefore becomes less prone to contractures.

Cycle duration: 3 weeks, 1 x/day

Program: **Capillarization 15** 

#### Muscular and circulation conditioning before physical activity

People who regularly practise sports are well aware of the transition, which is often painful, between a rest activity and an occasional physical effort.

The objective of the usual warming-up techniques is to respond to this need to activate progressively the physiological functions involved in one-time physical activity. This physiological need is also very desirable for moderate but unusual physical efforts, as is frequently the case for most of us (treks, biking, jogging, etc.).

The **Toning massage** program offers an ideal benefit in the form of optimal muscular and circulation conditioning before any type of physical effort. It allows us to avoid the uncomfortable sensations – usually experienced during the first few

minutes of unusual physical efforts – and to limit the secondary consequences of insufficient preparation (aches, etc.).

Cycle duration: To be used on the muscles most subject to stress in the physical activity involved (in this example: the quadriceps), in the last thirty minutes preceding physical activity

Program: **Toning massage** 8◇ 

**Use of Potentiation program to optimise the effects of explosive strength (sprints, jumps, squash, football, basketball, etc.) immediately before specific training or competition**

The **Potentiation** program should be used on the key muscles involved in the discipline practised. In this example, the key muscles of the sprinter (quadriceps) will be stimulated. For other disciplines the choice of muscles to be stimulated may be different.

The **Potentiation** program must not take the place of the voluntary warm-up usually performed before the competition. Activation of the cardiovascular system, short accelerations becoming progressively faster, practice starts and stretching will therefore be carried out by the athlete according to his normal routine.

A **Potentiation** session of short duration (approximately 3 minutes) will be applied on the sprinter's quadriceps immediately before the start of his race (or races, in the event of qualifying events). The specific muscular activation of the **Potentiation** program allows the maximum level of performance to be attained in the very first seconds of the race.

Cycle duration: 1 x

Always carry out the session as soon as possible before the start and at all events not more than 10 minutes before the start. After 10 minutes, the phenomenon of potentiation of muscle fibres quickly wears off

Program: **Potentiation** 8G 

Sport Category

Programs	Effects	Uses	Stimulation energies	Mi functions
<b>Endurance</b>	Improved absorption of oxygen by the stimulated muscles Improved performance for endurance sports	For athletes who wish to improve their performance in long-duration sporting events	Maximum tolerable energy (0-999)	Mi-SCAN
<b>Resistance</b>	Improved lactic capacity of muscles Improved performance for conscious resistance sports	For competitive athletes who wish to increase their capacity to obtain intense and prolonged effort	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
<b>Strength</b>	Increased maximum strength Increased rate of muscular contraction	For competitive athletes practising a discipline that requires strength and speed	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
<b>Explosive strength</b>	Increases the speed at which the level of strength is attained Improved efficiency of explosive actions (jumping, sprinting, etc.)	For athletes practising a discipline in which explosive strength is an important factor in performance	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
<b>Plyometry</b>	Improved strength/speed Work regime reproducing muscle activation induced by plyometric training	For athletes who, on account of injury, must interrupt or limit their voluntary plyometric training sessions  For athletes who wish to increase the quantity of plyometric muscular stimulations without risk of injury  To reproduce plyometric work on muscle groups for which this type of work is difficult to achieve voluntarily (deltoid, latissimus dorsi, etc.)	Maximum tolerable energy (0-999)  The energy is adjusted during phases of the program when the muscle reacts with muscle twitches	Mi-SCAN Mi-ACTION
<i>This program should be used under the supervision of a qualified sports trainer</i>				
<b>Hypertrophy</b>	Increased muscle volume Increased muscle resistance	For competitive athletes who wish to increase their muscle mass  For body-builders	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION

To define the level of Sport programs in relation with your personal characteristics, please see the training planner available on the DVD. At the end of a cycle, you may either start a new cycle at the next level up or do some maintenance at the rate of 1 session a week at the last level reached.

## Specific applications

### Preparation for a cyclist training three times a week who wishes to improve endurance performance

Effort expended over a long period of time makes demands on the aerobic metabolism, for which the decisive factor is the quantity of oxygen consumed by the muscles. To progress in endurance, it is therefore necessary to increase as much as possible the supply of oxygen to the muscles stimulated by this type of effort. Because oxygen is conveyed by the blood, it is essential to have an efficient cardiovascular system, due to voluntary training under certain conditions. However, a muscle capacity to consume the oxygen it receives (oxydative capacity) can also be improved by following a specific work regime.

The **Endurance** program of the Sport category leads to a significant improvement in the consumption of oxygen by muscles. Combining this program with the **Capillarization** program (of the Warm-up category), which develops the network of intramuscular capillaries, is particularly beneficial and allows endurance athletes to improve their performance levels.

Cycle duration: 8 weeks, 3 x/week

#### Eg for 1 week

Mon: Rest

Tue: 1 x **Endurance** 8G 

Wed: Cycle training 1h30 (moderate speed), then 1x **Capillarization** 8G 

Thu: 1 x **Endurance** 8G 

Fri: Rest

Sat: Cycle training 60' (moderate speed), then 1 x **Endurance** 8G 

Sun: Cycling 2h30 (moderate speed), then 1 x **Capillarization** 8G 

Programs: **Endurance** 8G and **Capillarization** 8G 

### Preparation for a runner training three times a week who wishes to progress in endurance (half-marathon, marathon)

Running as many miles as possible is essential in order to improve performance in endurance events. However, the strain on tendons and joints this type of training causes is today universally recognised. Integrating Complex electrostimulation into the training of the long-distance runner offers an excellent alternative that can help to overcome this problem. The **Endurance** program, which improves the muscles' ability to absorb oxygen, and the **Capillarization** program (of the Warm-up category), which develops the capillaries in the muscles, allow greater endurance to be achieved while limiting the weekly mileage and therefore the risk of injury.

Cycle duration: 8 weeks, 3 x/week

#### Eg for 1 week

Mon: Rest

Tue: 1 x **Endurance** 8G 

Wed: - Voluntary warm up training 20', then 1-2 series of 6 x (30" fast / 30" slow)

- Slow jogging 10' at the end of the session, then 1 x **Capillarization** 8G 

Thu: 1 x **Endurance** 8G 

Fri: Rest

Sat: Loose jogging 60', then 1 x **Endurance** 8G 

Sun: Extended run 1h30 (moderate speed), then 1 x **Capillarization** 8G 

Programs: **Endurance** 8G and **Capillarization** 8G 

### Pre-season preparation of lactic capacity for a resistance sport with three active training sessions per week (800 metres, track cycling, etc.)

Example of planning to develop the lactic capacity (resistance) of the quadriceps. For other disciplines, the choice of muscles to be stimulated may be different (to determine these muscles in relation to your sporting discipline, refer to the Training planner in the DVD).

During pre-season preparation for sports that make great demands on the lactic anaerobic system (intense effort sustained over as long as possible), it is essential not to neglect specific muscle preparation. Stimulation of the quadriceps (or a different key muscle group depending on the discipline practised) by means of the **Resistance** program results in improved anaerobic power, as well as greater muscle tolerance to high concentrations of lactates. The practical benefits are evident: improved performance thanks to better muscular resistance to fatigue for exercises of the lactic anaerobic type.

In order to optimise the effects of this preparation, you are advised to supplement it with **Capillarization** sessions, carried out during the week leading up to the competition (see "Warm-up category", "Achieving peak form before competition for resistance sports (eg: 800 metres, 1500 metres, etc.)").

An **Active recovery** (of the Recovery category) session is recommended after the most intensive training. It accelerates the rate of muscle recovery and reduces fatigue during the period of the season when the amount of training is high.

Cycle duration: 6-8 weeks, 3 x/week

#### Eg for 1 week

Mon: 1 x **Resistance** 8G 

Tue: Voluntary training on the track

Wed: 1 x **Resistance** 8G 

Thu: Voluntary intensive training on the track, then 1 x **Active recovery** 8G 

Fri: Rest

Sat: Rest

Sun: 1 x **Resistance** 8G , followed by session voluntary training on the track

Programs: **Resistance** 8G and **Active recovery** 8G 

### Preparation for a cyclist training three times a week who wishes to improve his power

Developing the strength of the thigh muscles is always beneficial for the competitive cyclist. Certain forms of training on the bicycle (hill work) can make a contribution in this respect. However, results will be more spectacular if additional training using Complex muscular stimulation is undertaken at the same time.

The special regime of muscle contractions of the **Strength** program and the large amount of work to which the muscles are subjected allow a significant increase of the strength of the thigh muscles.

Moreover, the **Active recovery** program (of the Recovery category), carried out within three hours of the most intensive training sessions, encourages muscular recuperation and makes it possible to follow on with qualitative training under optimum conditions.

Cycle duration: 8 weeks, 3 x/week

#### Eg for 1 week

Mon: Rest

Tue: 1 x **Strength** 8G 

Wed: - Cycle training 45' (moderate speed), then 5-10 times on a 500-700 m hill (rapidly)

- Recovery during descent

- Inactivity 15-20', then 1 x **Active recovery** 8G 

Thu: 1 x **Strength** 8G 

Fri: Rest

Sat: Cycle training 60' (moderate speed), then 1 x **Strength** 8G 

Sun: - Cycling 2h30-3 h (moderate speed)

- Muscle strengthening on hills (use of a high gear ratio remaining seated), then 1 x **Active recovery** 8G 

Programs: **Strength** 8G and **Active recovery** 8G 

### Preparation for a swimmer training three times a week who wishes to improve his swimming power

In swimming, developing the propulsive force of the upper limbs is an important factor in improving performance. Certain forms of voluntary training practised in the water can contribute to this. However, integrating Compex muscular stimulation into the voluntary training program makes it possible to achieve far better results. The special muscular contraction regime of the **Strength** program and the large amount of work to which muscles are subjected will allow you to increase significantly the strength of the latissimus dorsi, key muscles for the swimmer.

Moreover, the **Active recovery** program (of the Recovery category), carried out within three hours of the most intensive training, encourages muscular recuperation and makes it possible to follow on with qualitative training under optimum conditions.

Cycle duration: 8 weeks, 3 x/week

#### Eg for 1 week

Mon: Rest

Tue: 1 x **Strength** 18C 

Wed: - 20-30' swimming training (different styles), then 5-10 times 100 m with pull-boy

- Recovery 100 m backstroke

- Inactivity 15', then 1 x **Active recovery** 18D 

Thu: 1 x **Strength** 18C 

Fri: Rest

Sat: 1 h swimming training including some technical work, then 1 x **Strength** 18C 

Sun: - 20-30' swimming training (different styles), then 5-10 times 100 m with paddles

- Recovery 100 m backstroke

- Inactivity 15', then 1 x **Active recovery** 18D 

Programs: **Strength** 18C and **Active recovery** 18D 

### Pre-season preparation for a team sport (football, rugby, handball, volleyball, etc.)

Planning example to develop the strength of the quadriceps. Depending on the sport being practised, a different muscle group may be chosen.

During the preparatory pre-season period for team sports, it is essential not to neglect specific muscle preparation. In most team sports, the qualities of speed and strength make all the difference. Stimulation of the quadriceps (or a different key muscle depending on the sport concerned) by means of the Compex **Strength** program will result in increased speed of contraction and muscular strength. The practical benefits will be obvious: improved speed off the mark and during movement, jumping, shooting power, etc.

An **Active recovery** session (of the Recovery category), carried out after the most intensive training, accelerates muscular recuperation and reduces fatigue that has built up during the season when the workload is high.

Cycle duration: 6-8 weeks, 3 x/week

#### Eg for 1 week

Mon: 1 x **Strength** 8G 

Tue: Collective training, then 1 x **Active recovery** 8D 

Wed: 1 x **Strength** 8G 

Thu: Collective training, then 1 x **Active recovery** 8D 

Fri: 1 x **Strength** 8G 

Sat: Rest

Sun: Collective training or friendly match, then 1 x **Active recovery** 8D 

Programs: **Strength** 8G and **Active recovery** 8D 

### Maintaining results achieved during preparation for team sports during the competitive period (football, rugby, handball, volleyball, etc.)

This example concerns only athletes who have completed a full cycle of training by electrostimulation (at least 6 weeks) as part of their pre-season preparation. The weekly session of stimulation with the **Strength** program should be carried out on the same muscle groups as those stimulated during the preparation period (in our example, the quadriceps).

During the season, when matches are played regularly, care should be taken not to overtrain the specific musculature. On the other hand, nor should the benefits of preparation be lost by suspending stimulation training for too long a period. During this period of competition, muscle qualities should be maintained by means of one weekly session of stimulation carried out using the **Strength** program. It is also essential to leave a sufficiently long interval between this single weekly stimulation session and the day of the competition (3 days minimum).

The **Active recovery** program (of the Recovery category), which must be used during the first three hours after the match, as well as after each session of intensive training, restores muscle equilibrium more quickly.

Cycle duration: During the sport season, 1 x/week

#### Eg for 1 week

Mon: Rest

Tue: Collective training, then 1 x **Active recovery** 8D  (if training is intensive)

Wed: 1 x **Strength** 8G 

Thu: Collective training, then 1 x **Active recovery** 8D  (if training is intensive)

Fri: Rest

Sat: Rest

Sun: Match, then 1 x **Active recovery** 8D  (within the 3 hours that follow the competition)

Programs: **Strength** 8G and **Active recovery** 8D 

### Pre-season preparation of the explosive strength of the quadriceps for an athlete training three times a week (long jump or high jump, sprinting, etc.)

For other disciplines the choice of muscles to be stimulated may be different (refer if necessary to the Training planner in the DVD).

For all sports where the essential performance factor is explosive muscle strength, specific muscle preparation is the main factor in pre-season preparation. The explosive strength of muscles may be defined as the capacity of a muscle to attain a high level of maximum strength as quickly as possible. To develop this quality, voluntary training relies on tiring muscle training sessions that often include the risk of injury, since they are necessarily carried out using heavy weights. Integrating the use of the **Explosive strength** program lightens the muscle training sessions while at the same time offering greater benefits and more time for technical work.

Cycle duration: 6-8 weeks, 4 x/week

#### Eg for 1 week

Mon: 1 x **Explosive strength** 8G 

Tue: Voluntary training in stadium

Wed: 1 x **Explosive strength** 8G 

Thu: Voluntary training including technical work jumping pit

Fri: 1 x **Explosive strength** 8G 

Sat: Rest

Sun: Voluntary training in stadium followed by 1 x **Explosive strength** 8G

Program: **Explosive strength** 8G 

### Preparation for a thrower (javelin, shot putt) or tennis player (service, smash) who wishes to increase explosive strength through plyometric training

**Plyometric** training consists in combining an eccentric contraction with tensioning of elastic structures, followed by a rapid and explosive concentric contraction. During a sequence such as this, the rate of discharge of motor neurones is very distinctive and occurs with a very high frequency during the explosive concentric contraction.

The **Plyometry** program reproduces the sequence of nervous discharges during plyometric training, leading to development of the muscular structures as a result of this nervous stimulation. This high level of stimulation with the Compex increases the explosive strength of muscles, without the usual risk of injury connected with voluntary plyometric training.

Depending on the sport being practised, a different muscle group may be chosen (to determine this muscular group, refer to the training planner on the DVD). In this example, it will be beneficial for the javelin thrower and tennis player to stimulate the latissimus dorsi using the **Plyometry** program. The shot putt thrower or boxer, on the other hand, will work on the triceps.

Cycle duration: 9 weeks, 3 x/week

Progression in the levels:

Week 1: **Plyometry** level 1

Weeks 2-3: **Plyometry** level 2

Weeks 4-5: **Plyometry** level 3

Weeks 6-7: **Plyometry** level 4

Weeks 8-9: **Plyometry** level 5

Program: **Plyometry** 18C 

### Muscle volume for a body-builder

Despite repeated efforts during their voluntary training, many body-builders encounter difficulties in developing certain muscle groups. The specific stimulation imposed on muscles by the **Hypertrophy** program significantly increases the volume of the stimulated muscles. In addition, for a similar session time, the Compex **Hypertrophy** program provides a greater volume gain than voluntary training. The additional training imposed by this stimulation program on muscles not sufficiently receptive to traditional training provides a solution for the harmonious development of all muscle groups without recalcitrant areas.

To obtain optimum progress, you are advised to:

1) precede the **Hypertrophy** sessions with short voluntary training focused on strength; for example 3 series of 5 repetitions at 90% of maximum force;

2) carry out a **Capillarization** (of the Warm-up category) session directly after the **Hypertrophy** session.

### Muscle volume for a body-builder training three times a week

In this example, we assume that the body-builder wants to work on his/her biceps, which he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. biceps then calves).

Cycle duration: 8 weeks, 3 x/ week

Progression in the levels:

Week 1: **Hypertrophy** level 1

Weeks 2-3: **Hypertrophy** level 2

Weeks 4-5: **Hypertrophy** level 3


Weeks 6-8: **Hypertrophy** level 4

*Eg for 1 week*

Mon: Rest

Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 20D  followed by 1 x **Capillarization** 20D 

Wed: Rest

Thu: Voluntary training focused on the muscles of the trunk, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 20D  followed by 1 x **Capillarization** 20D 

Fri: Rest

Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 20D  followed by 1 x **Capillarization** 20D 

Sun: Rest

Programs: **Hypertrophy** 20D and **Capillarization** 20D 

### Muscle volume for a body-builder training five times a week

In this example, we assume that the body-builder wants to work on his/her calves that he/she considers to be his/her weak point. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. calves then biceps).

Cycle duration: 12 weeks, 5 x/ week

Progression in the levels:

Week 1: **Hypertrophy** level 1


Weeks 2-3: **Hypertrophy** level 2

Weeks 4-5: **Hypertrophy** level 3

Weeks 6-8: **Hypertrophy** level 4

Weeks 9-12: **Hypertrophy** level 5

*Eg for 1 week*



Mon: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then **Hypertrophy** 4A  followed by 1 x **Capillarization** 4D 

Tue: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 4A  followed by 1 x **Capillarization** 4D 

Wed: Rest

Thu: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 4A  followed by 1 x **Capillarization** 4D 

Fri: Voluntary training focused on the muscles of the lower limbs, followed by active work on the calf muscles: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 4A  followed by 1 x **Capillarization** 4D 

Sat: Voluntary training focused on the muscles of the upper limbs, followed by active work on the biceps: 3 series of 5 repetitions at 90% of Fmax., then 1 x **Hypertrophy** 4A  followed by 1 x **Capillarization** 4D 

Sun: Rest

Programs: **Hypertrophy** 4A and **Capillarization** 4D 

# Cross-training Category

Programs	Effects	Uses	Stimulation energies	mi functions
<b>Fartlek</b>	Training and preparation of muscles for all kinds of muscular work (endurance, resistance, strength, explosive strength) through different working sequences	At the beginning of the season to "re-set" the muscles after a rest period and before more intensive and specific training  During the season for those who do not have any preferences for a single kind of performance and prefer to exercise their muscles with different work-out methods	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
<b>Cross-endurance</b>	Training and preparation of muscles for all kinds of muscular work through different working sequences, concentrating on work of moderate intensity and long duration	For persons who want, above all, to improve their endurance, but also to develop other muscular qualities  Varied training concentrating on the development of aerobic muscular qualities	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
<b>Cross-resistance</b>	Training and preparation of muscles for all kinds of muscular work through different working sequences, concentrating on high-intensity work practised for as long as possible  Increased capacity of the fibres to support work leading to an accumulation of lactic acid	For persons who want to intensify their training schedule without increasing the risks of injury  To improve the resistance of muscular fibres to fatigue	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
<b>Cross-explosivity</b>	Training and preparation of muscles for all kinds of muscular work through different working sequences, with an emphasis on the development of contraction speed  Improvement of muscular explosive capacity	For persons who want to subject their muscles to different types of work, while concentrating on the improvement of contraction speed	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION

## Cross-training

### Specific applications

#### Resuming activity after a break in training: variable muscle preparation

In this example, we assume that the sportsman wants to work above all on the latissimus dorsi. It is of course possible to stimulate other muscles as well. Also this stimulation routine can be applied to several muscular groups at the same time (e.g. latissimus dorsi then quadriceps).

The fartlek comes from Scandinavia. It involves diversified training carried on in a natural setting. During training, different work-out routines are alternated and different muscles are stimulated. For example, after a few minutes of slow jogging, the trainee does a series of accelerations followed by some jumps, before resuming a slower jog, and so on. The objective is to work on different muscular qualities without however emphasising any one in particular. This activity is frequently done either at the beginning of the season, for general muscular reactivation, or done regularly by leisure-time sports enthusiasts who do not want to emphasise any particular kind of muscular performance, but rather want to stay fit and reach the right level in all kinds of muscular work.

At the beginning of the season or after a significant break in training, the resumption of physical and/or sports activity should be done progressively and become increasingly specific. It is thus normal to perform a few initial sessions with the aim of making the muscles do all kinds of work to prepare them for later training that will be more intensive and more oriented towards a specific kind of performance.

Through these eight sequences that automatically follow each other, the **Fartlek** program imposes different kinds of work to the stimulated muscles and thus makes them used to all kinds of effort.

Cycle duration: 1-2 week, 4-6 x/week

Program: **Fartlek** 18C 

#### Preparation for a long-distance runner who also wants to improve his speed at the end of the race

Any kind of physical activity of moderate intensity practised for a long period involves the aerobic metabolism. A good aerobic system depends on an efficient cardiovascular system. At the same time, it requires muscles capable of consuming the oxygen that reaches them (oxidative capacity). Electrostimulation develops this capacity and thus improves the endurance of the stimulated muscular groups.

However, it is useful, during a stimulation session, to change the focus of the program by alternating endurance with other aspects, such as speed, strength and resistance. This allows you not only to develop the endurance of your muscular fibres, but also to work simultaneously on other aspects that may be necessary when practising your chosen physical activity. A long-distance runner, for example, may need strength to cope with the slopes along the route or to increase his speed at the end of the race.

The characteristics of its sequences make the **Cross-endurance** program particularly suitable for improving endurance, but also for improving contraction speed, strength and resistance.

Cycle duration: 12 weeks, 3 x/week

Progression in the levels:

Week 1: **Cross-endurance** level 1

Weeks 2-3: **Cross-endurance** level 2

Weeks 4-6: **Cross-endurance** level 3

Weeks 7-8: **Cross-endurance** level 4

Weeks 9-12: **Cross-endurance** level 5

Eg for 1 week

Mon: Voluntary training: warm up 20', then 1-2 series of 6 x (30" fast / 30" slow)

Slow jogging 10' at the end of the session, then 1 session **Capillarization** 8 

Tue: 1 session **Cross-endurance** 8G 

Wed: Loose jogging 60', then 1 session **Capillarization** 8 

Thu: 1 session **Cross-endurance** 8G 

Fri: Rest

Sa : Loose jogging 45-60'

Di : Extended run 1h30 (moderate speed),  
then 1 session **Capillarization 8G** 

Programs: **Cross-endurance 8G** and  
**Capillarization 8G** 

### Improving the user's capacity for intense efforts without increasing the risks of injury

In this example, we assume that the person wants to work above all on his calves. It is of course possible to stimulate other muscles as well. Also, the **Cross-resistance** program can be used on several muscle groups at the same time, that is to say the same stimulation procedure can be applied to both calves and buttocks, for instance.

**Fitness enthusiasts** who regularly do several sessions of intensive voluntary activity every week frequently run the risk of over-stressing the muscles and joints involved. The **Cross-resistance** program offers a good alternative for intensifying the training schedule without increasing the risks of injuries. Also, the alternation of varied stimulation sequences means that contrasting types of work requiring several different muscular qualities can be imposed on the muscles, while also improving their capacity to support better the lactic acid accumulated during high-intensity physical effort.

Cycle duration: 3-6 sem., 3 x/sem.

#### Eg for 1 week

Mon: 1h of physical activity in the fitness centre

Tue: 1h of physical activity in the fitness centre, then 1 session **Cross-resistance 4A** 

Wed: Rest

Thu: 1h of physical activity in the fitness centre, then 1 session **Cross-resistance 4A** 

Fri: 1 session **Cross-resistance 4A** 

Sat: 1h-1h30 of voluntary physical activity (jogging, swimming, cycling)

Sun: Rest

## Cross-training

Program: **Cross-resistance 4A** 

### Competition preparation for an amateur practising a team sport

In this example, we assume that the person wants to work above all on his quadriceps. It is of course possible to stimulate other muscles as well. Also, the **Cross-explosivity** program can be used on several muscle groups at the same time, that is to say the same stimulation procedure can be applied to both quadriceps and buttocks, for instance.

**Team sports**, such as, for example, football, rugby, handball or basket-ball, require a good level of development of different muscular qualities. During a match, strength and explosive strength are decisive qualities, but a good cardio-vascular system plays an important part. Preparation for this type of physical activity should therefore not be concentrated exclusively on a single quality. Rather, it should take into account, through different exercises, the need for diversified voluntary preparation, which also involves factors such as muscular explosive strength, endurance and resistance.

**The Cross-explosivity** program activation routines allow varied work on the muscular groups involved and the selection of different working routines concentrating on muscular explosive strength, an important factor in most team sports.

Cycle duration: 10 weeks, 4 x/week

#### Progression in the levels:

Weeks 1-2: **Cross-explosivity** level 1

Weeks 3-4: **Cross-explosivity** level 2

Weeks 5-6: **Cross-explosivity** level 3

Weeks 7-8: **Cross-explosivity** level 4

Weeks 9-10: **Cross-explosivity** level 5

#### Eg for 1 week

Mon: Collective training, then 1 session **Cross-explosivity 8G** 

Tue: 1 session **Cross-explosivity 8G** 

Wed: Training 45'-1h in a fitness centre

Thu: Rest

## Cross-training

Fri: Collective training, then 1 session

**Cross-explosivity 8G** 

Sat: Rest

Sun: 1 session **Cross-explosivity 8G** 

Program: **Cross-explosivity 8G** 

# Conditioning Category

Programs	Effects	Uses	Stimulation energies	mi functions
<b>Abs</b>	Improvement of the maintenance of the abdominal wall muscles	During training or maintenance throughout the sports season	Maximum tolerable energy (0-999)	mi-SCAN
<b>Low back reinforcement</b>	Improvement of the active stability of the lower back region	During training or maintenance throughout the sports season	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
<b>Core stabilisation</b>	Improvement of the stabilisation of trunk muscles	During training or maintenance throughout the sports season	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION
<b>Ankle twist prevention</b>	Increase in the contraction speed and strength of the lateral peroneal muscles	To avoid recurrence of ankle sprain, for which its initial treatment permits restarting sports activities	Maximum tolerable energy (0-999)	mi-SCAN mi-ACTION

## Specific applications

### Abdominal muscle workout for golfers

Building the abdominal muscles is considered as a key in any general fitness program. Although rarely considered crucial to performance, the abdominal muscles have either primary or secondary importance for all sports activities. Many people performing abdominal exercises do not always understand the correct way to approach them and must be very careful not to over-curve the lower back (lordosis), which could lead to low back pain (lumbar pain). Isometric electrostimulation makes powerful muscle contractions possible without voluntary muscle movement, considerably reducing mechanical strain on the vertebrae and intervertebral discs of the lower back region. The **Abs** program makes it possible to perform a high quality workout perfectly suited to the characteristics of the abdominal wall muscles.

Cycle duration: 3-6 weeks, 3 x/week

Program: **Abs 49I-**

### Strengthening the lower back muscles in rowers

The lower back region consists of the final mobile vertebrae of the spinal column. Because of their physical position, these last vertebrae, as well as the surrounding structures (intervertebral discs, ligaments, and others), receive the greatest pressures, which explain the great frequency of lumbar pain, a health care phenomenon that gained particular prevalence during the 20th century.

The lower back muscles, along with the abdominal muscles, ensure the active maintenance of the lower part of the spinal column. When well toned, these muscles reduce the mechanical strain on the ligaments and discs of the lower back region, helping to prevent lumbar pain in concert with other preventive factors.

Certain sports activities, such as rowing, call upon significant use of the lower back region and, therefore, require specific training of the lower back muscles. These muscles act as a key link in the chain of muscles used for the specific motions of rowers and thereby contribute to the overall performance of the athlete.

The **Low back reinforcement** program specifically improves the quality of contractions of the lower back muscles. The program makes it possible to exercise these muscles in an isolated and intense manner to increase the active maintenance and strength of the lower back region.

Cycle duration: 3-6 weeks, 3 x/week

Program: **Low back reinforcement 12L-**

### Improving stabilisation of upper body muscles in triathletes

Most motions in sports activities involve one or more muscle chains under the control of a complex neuromuscular coordination and a specific technique, which both allow athletes to realize optimal mechanical output.

In this context, the stabilization of the trunk muscles makes it possible to achieve the most efficient transmission of forces exerted by the numerous muscles involved in any complex motion of any sports activity.

A physical activity, such as swimming, which involves a large number of lower and upper limb muscles, requires optimal muscle maintenance in order to obtain the highest possible level of performance.

Toning exercises focus on developing muscle aptitude to stabilise the lower part of the upper body and are part of the basic exercises in physical training; their goal is to improve the strength of the muscles of the abdomen: lower back and abdominal muscles.

The **Core stabilization** program makes it possible to use high quality electrostimulation to train these muscle groups during the same workout session.

Cycle duration: 3-6 weeks, 3 x/week

Program: **Core stabilization 12 and 34B-**

### Preventive stimulation of the lateral peroneal muscles to avoid or reduce the risk of recurrence of ankle sprain in footballers

Ankle sprain is one of the most prevalent sports injuries. It involves damage to the external ligaments of the ankle and occurs at different levels of seriousness. Most often, the injury is relatively benign, requiring, nonetheless, temporary interruption of the sports activity and, most importantly, predisposing the ankle to possible recurrence, the main complication, which can lead to the need for reconstructive surgery on the ligaments (ligament repair) in the long-term.

The ankle's anatomy exposes the joint to the risks incurred by forced foot movements that turn the foot inward (inversion). These risks become significant when this kind of movement, reaching beyond physiological limits, is realized at high speed, sometimes accompanied by a very large increase in strain resulting from the additional weight of an opponent, as it is often the case in contact sports, such as football.

Two small muscles, located on the external surface of the leg (small and short lateral peroneal muscles), ensure additional active protection under normal conditions by contracting with strength and speed as soon as there is a potentially harmful motion.

The main goal of traditional physical therapy for ankle sprains (after treatment for pain and inflammation) is to recover peroneal muscle competence, meaning the capacity to contract at maximum strength and speed when the ankle joint is subjected to a potentially harmful motion. This requires strong muscles, which the electrostimulation performed as part of the

**Ankle twist prevention** program helps to build in optimum fashion, and proprioception training that develops the muscle reflex that monitors the ankle joint.

The **Ankle twist prevention** program must only be used after initial treatment of pain and inflammation and only under medical supervision.

Electrostimulation sessions are carried out while standing. As a first step, the exercise is performed standing on both feet. Then, when it is possible, it is carried out standing only on the leg with the injured ankle, which is a difficult exercise that has the advantage of stimulating neuromuscular proprioception, the key goal of the exercise. By gradually advancing as part of a personalised program, it is possible to increase

## Conditioning

the level of difficulty of the electrostimulation session by adding the loss of visual control (closing the eyes) and changing the nature of the footing (working out on less and less stable supports).

Cycle duration: 2 to 4 weeks, 3 to 5 x/week

Week 5 and after: 1 x/every 10 to 14 days

Program: **Ankle twist prevention 2A-**

## Recovery Category

Programs	Effects	Uses	Stimulation energies	Mi functions
<b>Active recovery</b>	Marked increase in blood flow Accelerates the elimination of waste from the muscular contraction Endorphinic effect (see Pain category) Relaxing effect	To improve and accelerate the muscle recovery after an intensive exercise  To be used during the first 3 hours after each session of intensive training or after a competition	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Recovery plus</b>	Increases blood flow Fights pain Promotes muscle relaxation	To promote muscle recovery after intense physical activity that has led to temporary muscle failure  Performed one to three hours after the end of the physical activity	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Relaxing massage</b>	Decreased muscular tension Drainage of the toxins responsible for the exaggerated increase of muscular tone Effect of well-being and relaxation	To eliminate uncomfortable or painful sensations, following an exaggerated increase of muscular tone	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Reviving massage</b>	Considerable increase of blood flow in the stimulated region Improvement of tissue oxygenation Elimination of free radicals	To contrast effectively sensations of fatigue and localised heaviness	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Anti-stress massage</b>	Improvement of tissue vascularisation Decreased muscular tension associated with stress Return to calm	To restore a sensation of well-being after stressful situations	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Undulated massage 1</b>	Activation of the local circulation and improvement of exchanges at tissue level Reduction of unpleasant sensations	To eliminate uncomfortable body tensions and recover or achieve an optimal state of relaxation	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE

Programs	Effects	Uses	Stimulation energies	Mi functions
<b>Undulated message 2</b>	Improvement of venous return flow Relief of muscular tensions Sedative effect General relaxation	To recover sensations of well-being and relaxation	Increase the stimulation energy progressively until the most pleasant level is reached, while taking care, if possible, to obtain clearly visible muscular twitches or until a pronounced tickling sensation can be felt under the electrodes	Mi-SCAN Mi-RANGE
<b>Regeneration</b>	Analgesic effect through the release of endorphins Marked increase in blood flow encouraging oxygenation and drainage Activation of the oxydative metabolism Reactivation of proprioceptive pathways	To be used the day after competition as recovery training or as a complement to this type of training, the intensity of which can then be reduced	Increase the stimulation energies progressively until marked muscle twitches are obtained After 10 minutes, when tetanic contractions begin, the energies should be increased to the bearable maximum (0-999) After this phase of tetanic contractions, reduce the energy but ensure that muscle twitches continue to remain pronounced	Mi-SCAN Mi-ACTION
<b>Heavy legs</b>	Increased venous return Increased supply of oxygen to muscles Drop in muscular tension Elimination of tendency to have cramp	To remove the feeling of heaviness in the legs arising in unusual situations (standing for long periods, heat, hormonal imbalance linked to the menstrual cycle, etc.)	Increase the stimulation energies progressively until marked muscle twitches are obtained	Mi-SCAN Mi-RANGE
<b>Stretching</b>	To optimise the effects of the voluntary technique of stretching, by means of a reduction of muscular tonus obtained by specific activation of the antagonist of the stretched muscle (reciprocal inhibition reflex)	For all athletes who wish to maintain or improve their muscular elasticity To be used after all training sessions or during a specific stretching session	Increase progressively the intensities to obtain a tension that is clearly perceived	Mi-SCAN

Specific applications

Practical rules of use for the Undulated message 1 and Undulated message 2 programs

The Undulated message 1 and Undulated message 2 programs use a new stimulation mode that produces a particularly pleasant wave-like effect, thus increasing the effectiveness of these programs.

In most of the sequences of these two programs, stimulation is alternated on the four channels to reproduce faithfully different manual massage techniques (which must always be done in the direction of the venous return flow, that is to say from the extremities toward the heart). This is especially important to obtain the most benefit from the circulatory effects of these techniques and programs.

To take full advantage of this wave-like effect, particular care should be paid to the positioning of the electrodes, and the rules for their use indicated below should be followed.

**Rule n° 1:** number of stimulation channels to use

The correct use of the Undulated message 1 and Undulated message 2 programs requires the use of the 4 stimulation channels and, therefore, of 8 electrodes.

These 4 stimulation channels must be positioned:

- on a limb or a segment of a limb (e.g., a thigh),
- on the back,
- on two limbs (e.g., the two calves).

**Rule n° 2:** required order of the stimulation channels


The wave effect must necessarily follow the direction of the venous return flow. That is why the required order of the four stimulation channels must be strictly complied with.

The order of the stimulation channels is as follows:


- Channel 1: cable fitted with a transparent (Mi-SENSOR) or blue colour indicator.
- Channel 2: cable fitted with a green colour indicator,
- Channel 3: cable fitted with a yellow colour indicator,
- Channel 4: cable fitted with a red colour indicator.

It is always recommended to follow the general indications below:


Application on a limb or a segment of a limb

Application on the back

Application on the two limbs

### Increasing speed of recovery after intense competitive events, such as marathons

Derived from medical uses aimed at increasing blood flow, relieving pain, and reducing muscle stress, the **Compex Active recovery** program is probably the most widespread program used by athletes around the world.

Athletes practicing a wide variety of sports have recognised the effectiveness of this program.

Certain sports, such as marathons, bring about temporary muscle failure, most frequently exhibited by cramps that occur at the end of the effort.

Under those conditions, the **Recovery plus** program gives you results very similar to the classic **Active recovery** program, while reducing the risk of muscle cramps.

This must be in combination with other very simple measures, like drinking enough water, taking a shower, and delaying the use of electrostimulation until one or two hours after ending the activity.

Since the goal of the **Active recovery** and **Recovery plus** programs is not to develop muscle fitness, they can be used with less precise electrode placement than the classic physical training programs.

Under these conditions, the advantage is to be able to increase the stimulation of muscles or muscle groups without adversely affecting the desired results. An exhausted marathon runner, therefore, can benefit from the effectiveness of the technique for the quadriceps and calf muscles during a single session.

Cycle duration: 1 session one or two hours after the end of the activity that caused the muscle failure.

A second session about six hours after the end of the activity can increase the effectiveness of the technique.

Program: **Recovery plus 50 and 51** 

### Treatment of uncomfortable muscular tensions in the back of the neck

Remaining in a seated position for long periods, associated with repetitive movements of the upper limbs (as is frequently the case in front of a computer screen), is often responsible for an uncomfortable – or even painful – increase of muscle tension of the back of the neck muscles. Any other condition of stress may also create situations leading to a state of excessive muscular tension, which is often responsible for painful or uncomfortable sensations.

The in-depth effects produced by the **Relaxing massage** program enable to contrast effectively such painful sensations, with a particularly significant result in terms of relaxation.

Cycle duration: To be used on the back of the neck muscles, whenever an occasional sensation of painful muscular tension is felt; to be repeated, if necessary, in cases of particularly strong muscular tension

Program: **Relaxing massage 15** 

### Treatment of a localised sensation of heaviness or an occasional state of fatigue

The many stresses of everyday life are often responsible for uncomfortable, or even painful, physical sensations. A circulation slowdown is often the result of insufficient physical activity, frequently aggravated by the need in one's working life to remain in the same position for many hours (in a seated position, for example).

Although not serious, this simple “vascular slowdown” is, however, frequently the cause of unpleasant sensations (for example: a sensation of heaviness, often localised in the lower limbs, but also, at times, in another region of the body).

The **Reviving massage** program produces, in the greatest comfort, a reactivation of the blood circulation that allows an acceleration of tissue oxygenation and the elimination of painful sensations, due to insufficient physical effort.

Cycle duration: To be used on the calf muscles, whenever an occasional sensation of heaviness is felt; to be repeated, if necessary, in cases of persistent discomfort

Program: **Reviving massage 25** 

### Returning to a state of well-being after stressful situations

The many stressful situations of everyday life and the worries they entail create a certain degree of psychic tension to which each of us reacts differently. This hyperactivity of the brain and the emotions also has repercussions on our body, which may respond in many unwanted ways.

Undoubtedly, the most frequent problem is a sensation of discomfort associated with persistent muscular tension. This is sometimes intensified by physical constraints, such as the need to remain for a long time in positions unsuited to the well-being of our body. For a person under stress, it is difficult to relax. The muscles are always too tense, producing uncomfortable and even painful sensations.

In this kind of situation, the **Anti-stress massage** program can reduce muscular tensions through low-frequency stimulation producing a sensation of perfect calm.

Applied to the affected region, at the rate of two sessions a day, this program eliminates the unpleasant sensations produced by a stressful situation and produces a relaxing effect leading to a return to a state of well-being.

Cycle duration: Use this program daily during the stressful period or whenever a painful muscular sensation is felt

Program: **Anti-stress massage 15** 

### To eliminate unpleasant body sensations and recover an optimal state of relaxation

To relax and eliminate uncomfortable physical sensations are objectives pursued by many people, especially after a difficult and often tiresome day.

To achieve these objectives, the **Undulated massage 1** program offers an ideal solution for all those people who want to improve their well-being in this way. The automatic succession of the different program sequences allows the user to obtain beneficial selective effects to overcome the various habitual components of physical discomfort: circulatory slow-down, exaggerated muscular tone, micro-stresses on the joints, etc. Also, the wave-like effect obtained through the alternation of the stimulation phases on the 4 channels further increases the effectiveness of the program and produces a pleasant sensation that is particularly appreciated.

The program can be used effectively on all parts of the body; however, the back with the muscles that support the spinal vertebrae are particularly prone to these problems. The following indication presents a specific application to relieve

this region.

Cycle duration: Use this treatment at will, whenever you want to and/or need it

Program: **Undulated massage 1 45** 

### To restore the feeling of muscle equilibrium and soothe undesirable manifestations

Staying for a long time in uncomfortable static positions and/or continually performing tasks that involve a high level of physical stress, both associated with a state of psychological stress, are frequently responsible for physical "disorders".

Although generally speaking benign, all too frequently they entail psycho-bodily distress in our modern society.

Too much brain-work, usually referred to as "stress", and not enough or too much physical stress lead to this state of fatigue, tension or heaviness that produces a sensation of general discomfort.

The **Undulated massage 2** program offers a highly beneficial relaxation session for all those people who are faced occasionally or regularly with these negative situations. The judicious combination of the wave-like effect with the various stimulation sequences, each taking advantage of the recognized benefits of the latest advances in the medical field, allow the achievement of matchless results.

Cycle duration: Use this treatment at will, whenever you want to and/or need it

Program: **Undulated massage 2** 46 

### Use of the Regeneration program to eliminate muscle fatigue more quickly (cross-country running, football, basketball, tennis, etc.) and restore good muscular sensations more quickly

The **Regeneration** program, known also as the "day-after program" must be used on key muscles for the discipline being practised. In this example, key muscles for the cross-country runner (quadriceps) will be stimulated. For other disciplines, the choice of muscles to be stimulated may be different. Note that this type of program is particularly beneficial for all sports where competitions are repeated at frequent intervals, during tournaments and cups in various sports.

The stimulation session using the **Regeneration** program must be carried out the day after a competition to replace or supplement so-called "restoration" training, which can therefore be less intensive. Contrary to the **Active recovery** program, which provokes no tetanic contractions and which must be used during the three hours after the competition or intensive training, the **Regeneration** program is a form of light training which, in addition to an analgesic effect and an increased blood flow, aims to impose a small degree of anaerobic training and likewise to provoke slight tetanic contractions that are not tiring, making it possible to reactivate the proprioceptive pathways. Energy pathways are also gently stimulated, allowing their metabolic equilibrium to re-establish itself.

The session is composed of 6 stimulation sequences that follow on automatically:

- 1<sup>st</sup> sequence : analgesic effect
- 2<sup>nd</sup> sequence : marked increase in blood flow
- 3<sup>rd</sup> sequence : tetanic contractions to restore muscular sensation
- 4<sup>th</sup> sequence : activation of the oxydative metabolism
- 5<sup>th</sup> sequence : marked increase in blood flow
- 6<sup>th</sup> sequence : relaxing effect

Cycle duration: To be used throughout the season, according to the frequency of competitions; the day after each competition.

Program: **Regeneration** 8G 

### Prevention of the feeling of heaviness in the legs

The feeling of heaviness in the legs is due to a temporary deficiency in the venous return, and not to major organic lesions. Stagnation of the blood mass in the legs is encouraged in certain situations: prolonged periods in a standing position, long periods continuously seated, intense heat, etc. The insufficient supply of oxygen to the tissue, and in particular the muscles, that results from this situation, causes the feeling of heaviness and discomfort in the legs. Use of the **Heavy legs** program accelerates the return of blood to the veins and has an important relaxing effect on the painful muscles.

Cycle duration: Use this treatment whenever an occasional feeling of heaviness in the legs occurs

Program: **Heavy legs** 25 J-

### Preparation for a footballer wishing to optimise the effects of active stretching of the hamstrings muscles

The **Stretching** program consists in stimulating the antagonist (i.e. opposite) muscle of the muscle subjected to stretching in order to exploit a well-known physiological mechanism: reciprocal inhibition reflex. This reflex, based on the proprioceptive sensitivity of muscles, consists of a very marked muscular relaxation. This allows more efficient stretching, since it is carried out on a more relaxed muscle.

Stimulation is therefore performed on the muscle opposing the muscle being stretched. This stimulation consists of a progressive appearance and disappearance (slowly and lasting a long time) of contractions, with complete rest between contractions. It is during the contraction (which increases with each new level) that the athlete stretches the chosen muscle group by means of a traditional voluntary stretching technique.

In this example, the stimulation is carried out on the quadriceps in order to facilitate stretching of the hamstrings muscles during the contraction phases.

Cycle duration: Throughout the season, according to the frequency of voluntary stretching sessions

According to the normal duration of stretching (this depends on schools of physical training and each person's feelings), choose the level that seems most appropriate. The proposed duration of stretching is:

- 10 seconds for level 1
- 12 seconds for level 2
- 14 seconds for level 3
- 16 seconds for level 4
- 18 seconds for level 5

Program: **Stretching** 8 

Adopt the starting position for the active stretching exercise

# Pain relief Category

Table of the pathologies

Pathologies	Programs	References
<b>Neuralgia of the upper limb</b> (brachial neuralgia)	<b>Modulated TENS</b>	See application page 148
<b>Chronic muscular pain</b> (polymyalgia)	<b>Endorphinic</b>	See application page 148
<b>Contracture</b> (eg: localized contracture in external side of the calf)	<b>Decontracture</b>	See application page 148
<b>Chronic muscular pain in the back of the neck</b> (cervical pain)	<b>Cervical pain</b>	See application page 148
<b>Muscular pain in the thoracic region</b> (thoracic back pain)	<b>Thoracic back pain</b>	See application page 149
<b>Muscular pain in the low back region</b> (low back pain)	<b>Low back pain</b>	See application page 149
<b>Sharp and recent muscular pains affecting a muscle in the low back</b> (lumbago)	<b>Lumbago</b>	See application page 149
<b>Chronic elbow pain</b> (epicondylitis = tennis elbow)	<b>Epicondylitis</b>	See application page 149
<b>Recent and acute back of the neck muscle pains</b> (torticollis)	<b>Torticollis</b>	See application page 150

## Pain relief

Programs	Effects	Uses	Stimulation energies	mi functions
<b>Modulated TENS</b>	Blocks transmission of pain by the nervous system	All acute or chronic localized pain	Progressively increase the stimulation energies until a pronounced tickling sensation can be felt under the electrodes	mi-TENS
<b>Endorphinic</b>	Analgesic action through the release of endorphins Increased blood flow	To combat chronic muscular pain	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Decontracture</b>	Reduced muscular tension Relaxing effect	To combat recent and localized muscular pain	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Cervical pain</b>	Analgesic action through the release of endorphins Increased arterial flow	Analgesic current specifically adapted to pains in the back of the neck	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Thoracic back pain</b>	Analgesic action through the release of endorphins Increased arterial flow	Analgesic current specifically adapted to pain in the thoracic region (below the back of the neck and above the low back)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Low back pain</b>	Analgesic action through the release of endorphins Increased blood flow	Analgesic current specifically adapted to persistent low back pain (lumbar region)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Lumbago</b>	Reduced muscular tension Relaxing effect	Analgesic current specifically adapted to sharp and sudden low back pain (lumbar region)	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE
<b>Epicondylitis</b>	Blocks transmission of pain by the nervous system	Analgesic current specifically adapted to persistent pain in the elbow	Progressively increase the stimulation energies until a pronounced tickling sensation can be felt under the electrodes	mi-TENS
<b>Torticollis</b>	Reduced muscular tension Relaxing effect	Analgesic current specifically adapted to acute and sharp pains in the back of the neck	Progressively increase the stimulation energies until pronounced muscle twitches are obtained	mi-SCAN mi-RANGE

■ The use of Pain relief programs should not be prolonged without medical advice. If the pain is great and/or persistent, you are advised to consult a doctor. Only a doctor can make an accurate diagnosis and put in place therapeutic actions designed to encourage the disappearance of the disorder.

## Specific applications

### Neuralgia of the upper limb (brachial neuralgia)

Some people suffer from arthritis in the joints of the vertebrae at the back of the neck, or from arthritis or periarthritis of the shoulder. These situations often give rise to pain that runs down one arm and is known as "brachial neuralgia". These arm pains that start in the shoulder or the back of the neck can be reduced with the Compex **Modulated TENS** program by following the practical recommendations outlined below.

Cycle duration: 1 week, 1 x/day minimum, then adapt according to how the pain develops

According to requirements, the Modulated TENS program can be repeated a number of times during the same day

Program: **Modulated TENS** 35◊

### Chronic muscular pain (polymyalgia)

Some people suffer from muscular pain that affects several muscles or parts of different muscles at the same time. The localization of this chronic pain can vary over time. These continual and diffuse muscle pains are the result of chronic contractures in which acids and toxins accumulate, irritating the nerves and causing pain. The **Endorphinic** program is particularly effective against such pain since, in addition to its pain-reducing effect, it increases the blood flow to the contracted muscular groups and removes from them accumulations of acids and toxins.

**Example:** localized pain in the biceps. **However,** the pain may affect other muscular groups. The practical application illustrated below remains valid, but it is then necessary to place the electrodes on the muscular group concerned.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Endorphinic** 20◊ 

### Contracture (eg: localized contracture in external side of the calf)

After tiring muscle work, intense training or a sporting competition, certain muscles or certain parts of muscles often remain tense and slightly painful. These are known as muscular contractures which should disappear after a few days with rest, good rehydration, a balanced food intake with mineral salts and application of the **Decontracture** program. The phenomenon of contracture frequently affects the calf muscles, but can also occur with other muscles. In this case, simply follow the same practical recommendations as below, but place the electrodes on the muscular group concerned.

Cycle duration: 1 week, 1 x/day

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Decontracture** 24◊ 

### Chronic muscular pain in the back of the neck (cervical pain)

Positions in which the muscles of the back of the neck remain tense for long periods of time, for example working in front of a computer screen, may be responsible for the onset of pain in the back of the neck or on either side of the base of the neck, in the upper back. These pains are due to a contracture of the muscles, of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If this phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or appears after only a few minutes spent working in an unfavourable position.

**These chronic pains in the back of the neck can be treated effectively with the Cervical pain program,** which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Cervical pain** 15◊-

### Muscular pain in the thoracic region (thoracic back pain)

Vertebral arthritis and positions in which the muscles of the spinal column remain under tension for long periods of time are often responsible for the onset of pain in the middle of the back, which is accentuated with fatigue. Pressure of the fingers on the muscles on either side of the spinal column can often trigger a sharp pain. **These pains are due to a contracture of the muscles,** of which the prolonged tension crushes the blood vessels and prevents the blood from supplying and oxygenating the muscle fibres. If the phenomenon is prolonged, there is an accumulation of acid and the blood vessels atrophy. The pain then becomes continuous or can appear after only a few minutes spent working in an unfavourable position.

**These chronic back pains can be effectively treated with the Thoracic back pain program** which reactivates the circulation, drains accumulations of acid, oxygenates the muscles, develops the capillaries and relaxes the contracted muscles.

Cycle duration: 4 weeks, 2 x/week with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Thoracic back pain** 13◊-

### Muscular pain in the low back region (low back pain)

Low back pain is the most frequently encountered pain. In a standing position, the entire weight of the trunk is concentrated on the joints between the last vertebrae and the sacrum. The low back region is therefore under particular strain. The discs between the vertebrae are crushed and the low back muscles contracted and painful. **There are very many treatments for the relief of low back pain sufferers;** among these, the specific currents of the Compex **Low back pain** program provide an appreciable improvement and can even

resolve the problem if it is essentially muscular in origin.

Cycle duration: 4 weeks, 2 x/day, with a 10 minute break between the 2 sessions

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Low back pain** 12◊-

### Sharp and recent muscular pains affecting a muscle in the low back (lumbago)

**During back movement,** for example when lifting something, when turning or when standing up straight after bending down, a sudden pain can be triggered in the low back. Those who suffer from this problem present a contracture of the low back muscles and feel a sharp pain in this region; because they cannot stand completely straight, they remain bent over on one side. All of these symptoms indicate what is known as lumbago, which is principally the result of a sharp and intense contracture of the low back muscles (lumbar region).

**In a situation such as this,** it is always necessary to consult a doctor to receive appropriate treatment.

**In addition to such treatment,** the specific Compex **Lumbago** program can help effectively to relax the muscles and remove the pain.

Cycle duration: 1 week, 1 x/day

You are advised to consult your doctor if no improvement is observed after the first week of use

Program: **Lumbago** 33◊-

### Chronic elbow pain (epicondylitis = tennis elbow)

All tendons of the muscles that allow us to stretch our hand, wrist and fingers end at the small external bone mass of the elbow (epicondyle). Hand and finger movements therefore transmit tensions that are concentrated in the tendon endings at this bone mass.

When hand movements are repetitive, as is the case for painters, tennis players or even those who constantly use the mouse of a computer system, small injuries, accompanied by inflammation and pain, develop in the region of the epicondyle. This is known as "epicondylitis", which is characterised by pain around the external bone mass of the elbow when pressure is applied or when the forearm muscles are contracted.

The **Compex Epicondylitis** program supplies specific current to combat this type of pain. It acts effectively as a supplement to rest. However, it is necessary to consult your doctor if the pain gets worse or does not disappear quickly after a few sessions.

Cycle duration: 1 week, 2 x/week minimum, then adapt according to how the pain develops

According to requirements, the Epicondylitis program can be repeated several times during the same day

Program: **Epicondylitis** 36↵

## Pain relief

### Torticollis (recent and acute back of the neck muscle pains)

As a consequence of various circumstances, such as, for example, a rather sudden movement of the neck, an incorrect sleeping position or exposure to a cold environment, a sharp pain may be felt in the back of the neck often accompanied by a considerable reduction of the mobility of the cervical region. A person affected by this problem has contracted muscles in the back of the neck and feels a sharp pain there. Together these symptoms represent what is known as "torticollis", which is mainly the result of a sharp and intense contracture of the muscles in the back of the neck.

In this kind of situation, it is always necessary to consult a doctor in order to receive appropriate treatment.

In addition, **Compex'** specific **Torticollis** program can effectively help to relax the muscles and eliminate the pain.

Cycle duration: 1 week, 1 x/day

You are advised again to consult your doctor if no improvement is observed after the first week of use

Program: **Torticollis** 15↵

## Rehabilitation Category

Programs	Effects	Uses	Stimulation energies	Mi functions
<b>Disuse atrophy</b>	Reactivation of the trophic action of muscle fibres damaged during atrophy	Treatment to combat any reduction of muscle volume: - as a result of trauma necessitating immobilisation - accompanying degenerative lesions to joints	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
<b>Muscle growth</b>	Increased diameter and capacity of muscle fibres damaged during under-activity or inactivity imposed by some kind of pathological problem Restoration of muscle volume	After a period of using the <b>Disuse atrophy</b> treatment, as soon as the muscle shows signs of a slight recovery of volume or tone Until the virtually complete restoration of muscle volume	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION
<b>Reinforcement</b>	Increase the strength of a muscle previously atrophied Increase the strength of a muscle affected by a pathological process	At the end of rehabilitation, when the muscle has regained normal muscle volume From the onset of rehabilitation for non-atrophied muscles	Maximum tolerable energy (0-999)	Mi-SCAN Mi-ACTION

In a rehabilitation process, it is essential, initially, to treat the atrophy (i.e. regain normal muscle volume) before attempting to increase the strength of the muscle using the **Reinforcement** program.

The use of programs of the **Rehabilitation** category must under no circumstances replace rehabilitation sessions carried out in the presence of the physiotherapist. Although the restoration of initial muscular qualities is a key factor in the process of rehabilitation, other aspects (joint mobility, vigilance, residual pain, etc.) can only be dealt with effectively by a competent health professional.

Some pathologies and some post-operative rehabilitations require special precautions during use; you are therefore always advised to seek the opinion of your doctor or physiotherapist before using any of the programs of the **Rehabilitation** category.

■ At the end of a cycle, you may do some maintenance at the rate of 1 session a week at the last level reached.

## Specific applications

### Atrophy of the quadriceps as a result of trauma

The quadriceps is a voluminous muscle located in the anterior part of the thigh. It is the main muscle that allows you to stretch the knee; its role is therefore essential for walking, running, climbing stairs, etc. Any trauma affecting a lower limb therefore results in a wasting of this muscle, the reduction of muscle volume being more or less serious depending on the duration of the period of inactivity.

This atrophy is normally spectacular when it occurs after trauma to the knee, particularly if the trauma was treated by means of surgical intervention.

The programs of the **Rehabilitation** category are specifically designed to treat the deterioration of muscle fibres that results from such a process. The progressiveness of the work imposed by the different recommended programs is decisive in obtaining optimum results.

Cycle duration: 10 weeks, 1 x/day

Weeks 1-2: **Disuse atrophy**

Weeks 3-8: **Muscle growth**

Weeks 9-10: **Reinforcement**

Programs: **Disuse atrophy**, **Muscle growth** and **Reinforcement** 8G 

### Development of the latissimus dorsi to treat and prevent tendinosis shoulder pain (rotator cuff syndromes)

The shoulder is a complex joint enabling us to make broad gestures (for example, raising our arms in the air). During some of these movements, the tendons in the shoulder can rub against or be compressed against bony fragments of the joint. When this phenomenon is repeated, or occurs in certain patients who have an unfavourable anatomic constitution, such damage to the tendons causes them to become inflamed and thicker, resulting in a considerable increase in their degree of compression. The pain often then becomes very severe, preventing all movement of the shoulder; it can even set in at night and cause serious sleep deprivation.

Appropriate medical treatment can only be put in place by consulting a doctor.

However, electrostimulation of the latissimus dorsi by means of specific Compex programs can reduce the distress to which tendons are subjected by increasing their freedom of movement around the shoulder joint.

Cycle duration: 6 weeks, 1 x/day

Weeks 1-2: **Disuse atrophy**

Weeks 3-6: **Muscle growth**

Programs: **Disuse atrophy** and **Muscle growth** 18C 

### Development of the abdominal belt to prevent pain in the lumbar region (low back pain)

Low back pain arises most frequently among subjects who do not have sufficient musculature in the abdominal region.

These muscles in fact represent a veritable natural "corset", the role of which is to protect the low back region from excessive stress of all kinds. This is why it is referred to as the abdominal "belt". After a lumbago, when the pain has stopped, a common recommendation is to improve the efficiency of the abdominal and back muscles to prevent any recurrence. The **Disuse atrophy** program imposes a large amount of work on the abdominal muscles, without requiring the harmful or even dangerous positions often adopted when voluntary exercises are carried out incorrectly. An abdominal belt that has greater strength and endurance can then satisfactorily fulfil its role of protecting the low back region.

Cycle duration: 4 weeks, 1 x/day

Program: **Disuse atrophy** 10I 

### Development of the low back muscles to prevent pain in the lumbar region (low back pain)

Like the muscles of the abdominal region, the muscles of the low back (lumbar muscles) also play a role in protecting the low back region. Subjects whose low back muscles are inadequately effective are particularly prone to low back pain. Once the pain in the low back has disappeared, patients are often advised to strengthen

the low back muscles with a view to preventing the recurrence of painful episodes. However, carrying out voluntary back exercises often presents serious difficulties for patients who suffer from low back pain. That is why electrostimulation of the low back muscles using the **Disuse atrophy** program is one of the preferred methods to improve the efficiency of these muscles.

Cycle duration: 4 weeks, 1 x/day

Program: **Disuse atrophy** 14B 

### Strengthening of the lateral peroneus muscles after ankle sprain

The purpose of the lateral peroneus muscles is to maintain the stability of the ankle joint and prevent it from rotating inwardly. After a sprain, these muscles lose their reflex-contraction capacity together with much of their strength. Regaining competent lateral peroneal muscles after a sprain is a fundamental step, without which recurrence is very probable. To do their job correctly, the lateral peroneals must be strong enough to prevent the foot twisting inwards, but they must also contract reflexively at the precise moment when the heel tilts inwards. To develop both of these aspects, strength and speed of contraction, you should use the **Reinforcement** program, which produces efficient lateral peroneal muscles and therefore helps to prevent recurrence.

Cycle duration: 4 weeks, 1 x/day

Weeks 1-4 : **Reinforcement**

Weeks 5 and following weeks: **Reinforcement** 1x/week if you take part in a dangerous sporting activity

Program: **Reinforcement** 2A 

# Test Category

Programs	Effects	Uses	Stimulation energies	mi functions
<b>Motor point</b>	Stimulation at low frequencies allowing the extremely precise identification of the ideal position for placing the positive electrode (motor point)  The precise positioning of the positive electrode on the motor point of the muscle to be stimulated is an important factor for the effectiveness of treatment	Recommended before any initial muscular stimulation session  The motor points thus precisely identified may be easily located again if they are marked using a skin marker pencil	Progressively increase the energy, while manually moving the pen tip on the conductive gel, without ever interrupting skin contact  A maximum energy between 5 and 15 is usually sufficient to identify the motor point precisely	None
<b>Demo</b>	Refer to the quick start guide, "Test Your Compex in 5 Minutes"			mi-SCAN mi-TENS mi-ACTION mi-RANGE

## Specific applications

### Locating a motor point (leg: locating the motor points of the vastus internus and the vastus externus of the quadriceps)

Electrostimulation subjects the muscles to work. The progress achieved depends on the kind of work to which the muscles are subjected, that is to say the program chosen. The electrical pulses generated by these programs are transmitted to the muscles (via the motor nerve) through self-adhesive electrodes. The positioning of the electrodes is one of the determinant factors for ensuring a comfortable electrostimulation session.

It is therefore essential to dedicate special care to this aspect. The correct placement of the electrodes and the use of significant energy allow a large number of muscle fibres to work. The greater the energy, the greater the spatial involvement, that is to say the number of fibres working, and therefore the greater the number of fibres that make progress.

A stimulation lead wire consists of two electrodes:

- a positive electrode (+): red connection,
- a negative electrode (-): black connection.

The positive electrode must be fixed precisely on the motor point of the muscle. The motor points correspond to an extremely localised area where the motor nerve is more excitable. Although the location of the various motor points is now well known, there may nevertheless be variations, which can extend to several centimetres, between different individuals.


The **Motor point** program, associated with the use of the motor point pen, allows to determine with great accuracy the exact location of the motor points for each individual and thus ensure the greatest effectiveness of the programs.

It is advisable to use this program before any initial muscular electrostimulation session. The motor points once located can be easily identified using a skin-marker pencil or in any other way, thus avoiding the need to repeat this search before each session.

## Advice on use

To locate the motor point of the vastus internus, proceed as described below:

- 1) Apply a large electrode at the top of the thigh (see picture n° 7),
- 2) Connect the negative terminal (black connection) of channel 1 to the output of this large electrode located towards the internal surface of the thigh,
- 3) Spread a thin but uniform layer of Compex conductive gel in the position indicated for the positive electrode positioned over the vastus internus (see picture n° 7), spreading the gel a few extra centimetres in all directions,
- 4) Connect the positive connection (red) to the tip of the motor point pen and bring the tip of the pen in contact with the conductive gel,

- 5) Select the **Motor point** program and set the area for the thighs; then, start the program by pressing the "+" / "-" keys of channel 1,


6) Very progressively increase the energy of channel 1, until a value between 5 and 15 is reached, while continuously moving the pen tip over the gel layer, but without ever losing contact with the skin gel, to avoid triggering off an electrode fault message,

7) As soon as you observe a muscular response in form of twitching, you have located the vastus internus motor point. Visually locate this motor point and apply a small electrode that should be centred over the motor point. If this causes an electrode fault, ignore the message and do not exit the program: continue with the rest of the operations normally,

8) Remove the pen from the positive connection (red) of channel 1 and connect the small electrode to this connection.

To locate the motor point of the vastus externus, proceed as described below:

- 1) Connect the negative terminal (black connection) of channel 2 to the other output of the large electrode that should be left in place,
- 2) Spread a thin but uniform layer of Compex conductive gel in the position indicated for the positive electrode positioned over the vastus externus (see picture n° 7), spreading the gel a few extra centimetres in all directions,
- 3) Connect the positive connection (red) to the tip of the motor point pen and bring the tip of the pen in contact with the conductive,

- 4) Select the **Motor point** program and set the area for the thighs; then, start the program by pressing the "+" / "-" keys of channel 2,

5) Very progressively increase the energy of channel 2, until a value between 5 and 15 is reached, while continuously moving the pen tip over the gel layer, but without ever losing contact with the skin gel, to avoid triggering off an electrode fault message,

6) As soon as you observe a muscular response in form of twitching, you have located the vastus externus motor point. Visually locate this motor point and apply a small electrode that should be centred over the motor point. If this causes an electrode fault, you can exit the program.

7) Remove the pen from the positive connection (red) of channel 2 and connect the small electrode to this connection.

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